



CALIFORNIA DEVICE RECALL INFORMATION SHEET

Hamilton Medical Recall Hamilton C-6 Ventilators For Sensor Problems

Recall Date	Product Description	Recalling Firm	Recall Reason
06/26/2024	Ventilator HAMILTON-C6 PN: 160021	HAMILTON MEDICAL AG Bonaduz, Switzerland	Ventilator may enter sensor fail mode, ventilation may not re-initiate, after patient is reconnected, which may cause hypoxia, if this sequence occurs: 1)User presses O2 enrichment, 2)User disconnects endotracheal tube for open suctioning, 3)Sensor error initiated, sensor fail mode occurs, 4)Patient is reconnected with sensor fail mode active, 5) Ventilation is not re-initiated by the ventilator.

Recall Class	Product Identification	Distribution	Affected Dates
I	Ventilator Hamilton-C6 UDI: 07630002808590, Software Versions: SW v1.1.4, SW v1.1.5, SW v1.1.6	68 Units Nationwide	May 2024 and prior

FOR ADDITIONAL INFORMATION, PLEASE VISIT THE [FDA WEBSITE](#)

CDPH Food and Drug Branch
MS 7602 • P.O. Box 997435 • Sacramento, CA 95899-7435
(916) 650-6500 • (916) 650-6650 FAX
Internet Address: www.cdph.ca.gov

