



**CALIFORNIA DEVICE RECALL INFORMATION SHEET**

**Philips North America Recalls Sense XL Torso Coil for Intera and Achieva MR**

Recall Date	Product Description	Recalling Firm	Recall Reason
07/03/2024	<b>Sense XL Torso Coil 1.5T</b> 16-element receive only coil for torso and abdomen imaging in Philips MR Intera and Achieva Systems	<b>PHILIPS NORTH AMERICA LLC</b> Cambridge, Massachusetts	Potential for coils to heat up and harm patients (burn).
07/03/2024	<b>Sense XI Torso Coil 3.0T</b> 16-element receive only coil for torso and abdomen imaging in Philips MR Intera and Achieva Systems	<b>PHILIPS NORTH AMERICA LLC</b> Cambridge, Massachusetts	Potential for coils to heat up and harm patients (burn).
07/03/2024	<b>Sense XL Torso Coil 1.5T mk 2</b> 16-element receive only coil for torso and abdomen imaging in Philips MR Intera and Achieva Systems	<b>PHILIPS NORTH AMERICA LLC</b> Cambridge, Massachusetts	Potential for coils to heat up and harm patients (burn).

Recall Class	Product Identification	Distribution	Affected Dates
I	<b>Sense XL Torso Coil 1.5T</b> Model No. 453567141882, 453567141883, 989603014351 & 989603014352; UDI-DI: N/A; Serial No. 288 ... 1314; 300202352, 336 ... 1319	729 Units Nationwide	June 2024 and prior
I	<b>Sense XI Torso Coil 3.0</b> Model No. 453567394942, 453567394945, 453567394943, 989603050642 & 989603050641; UDI-DI: N/A; Serial No. 529 ... 341	307 Units Nationwide	June 2024 and prior

I	<b>Sense XL Torso Coil 1.5T mk 2</b> Model No. 453567394942, 453567394945, 453567394943, 989603050642 & 989603050641; UDI-DI: N/A; Serial No. 16 ... 40000	14 Units Nationwide	June 2024 and prior

FOR ADDITIONAL INFORMATION, PLEASE VISIT THE [FDA WEBSITE FOR 1.5T](#) , [3.0T](#) and [1.5T MK2](#).